

Gel Extraction (Qiagen)

Written by
Tuesday, 30 November 2010 10:53 -

- Excise the DNA fragment from the Agarose gel with a clean, sharp scalpel.

(Minimize the size of the gel by removing extra agarose).

- Weigh the gel slice in a colourless tube. Add 3 volumes of Buffer QG to 1 volume of gel (100 mg~300µl).

- Incubate at 50°C for 10 mins. to dissolve agarose. Vortexing in between helps (37°C is fine).

- Add 1 gel volume of iso-propanol to the sample and mix (not for 500bp-4kb).

- Place a QIAquick spin column in a provided 2 ml collection tube.

Gel Extraction (Qiagen)

Written by

Tuesday, 30 November 2010 10:53 -

- To bind DNA, apply the sample to the QIAquick column and centrifuge for 1 min (maximum volume is 800µl, so do twice).

- Discard flow-through and place QIAquick column back in the same collection tube.

- Add 750µl Buffer PE to the column and centrifuge for 1 min.

- Discard flow-through and centrifuge the column for an additional 1 min. at 13,000 rpm.

- Place the column into a clean 1.5ml eppendorf.

Gel Extraction (Qiagen)

Written by

Tuesday, 30 November 2010 10:53 -

- To elute DNA, add 30µl Buffer EB or H₂O to the centre of the QIA membrane and centrifuge for 1 min.

- Can store at -20°C.